



**NICHIRYO** THE PIONEER IN FLUID SAMPLING AND DISPENSING SYSTEMS

# ACCUPENSER™

## JUNIOR

1ml 2ml 5ml 10ml 30ml

### THE MOST ACCURATE DISPENSER

- Dispensing mechanism is chemically inert, permitting use with all liquids except hydrofluoric acid and nitric acid
- Autoclavable without disassembling
- Graduated scale can be moved to lie horizontally without altering volume setting. This reduces the height of the dispenser which facilitates storage in a refrigerator
- Priming without loss of reagent

