



NICHIRYO THE PIONEER IN FLUID SAMPLING AND DISPENSING SYSTEMS

MODEL 2100DG

1 ml 2 ml 5 ml 10 ml

accuracy < $\pm 1.0\%$
precision < $\pm 0.2\%$

A DIGITALLY VOLUME READING HANDY DISPENSER

- EXCELLENT ACCURACY AND REPRODUCIBILITY
- DISPENSING FROM MOST LIQUID CONTAINERS IS AVAILABLE.
- VOLUME ADJUSTMENT IS SIMPLE AND REMAINS FIXED.
- LIGHT WEIGHT DESIGN PERMITS ONE HAND OPERATION.
- ALL PARTS WHICH MAY CONTACT DIRECTLY WITH LIQUID ARE AUTOCLAVABLE.

